

# ONE 3D Printing Solution

## Technical Specifications

| ONE 3D Printer                                  |                                                                             |
|-------------------------------------------------|-----------------------------------------------------------------------------|
| Technology                                      | 385nm DLP with Force Feedback                                               |
| Build Plate Size (W x D x H)                    | ~ 2 die models flat<br>5.2 x 3 x 4 inches                                   |
| Rapid Print Kit<br>Build Plate Size (W x D x H) | Up to 2 4-unit bridges<br>1.6 x 0.8 x 4 inches                              |
| XY Accuracy                                     | +/- 34µm                                                                    |
| Material library                                | Selected material availability                                              |
| Touch Screen Size                               | 7 inch touch screen                                                         |
| Connectivity                                    | WiFi, Ethernet, USB                                                         |
| Device Dimension                                | 12.6 x 14.8 x 15.4 inches                                                   |
| Device Weight                                   | 8 kg / 17.6 lbs.                                                            |
| Warranty                                        | 1-year CARE Plus Plan with remote onboarding and continued support included |



| WASH Post-Processing       |                                              |
|----------------------------|----------------------------------------------|
| Cleaning Mechanism         | Pump system with 2x propellers               |
| Chamber Volume (W x D x H) | max. 6 dental models<br>4.7 x 6.3 x 4 inches |
| Cleaning Agent             | Isopropylalcohol, Ethanol                    |
| Cycle Time                 | 6 - 10 minutes per cycle                     |
| Connectivity               | WiFi, Ethernet                               |
| Device Dimension           | 13.4 x 18.2 x 11.6 inches                    |
| Device Weight              | 15 kg / 33 lbs.                              |
| Warranty                   | 1-year CARE Plus Plan included               |

| CURE Post-Processing       |                                            |
|----------------------------|--------------------------------------------|
| Curing System              | 360° LEDs                                  |
| Chamber Volume (W x D x H) | max. 7 dental models<br>7 x 4.6 x 4 inches |
| Curing Environment         | Vacuum or inert gas                        |
| Cycle Time                 | 6 - 10 minutes per cycle                   |
| Connectivity               | WiFi, Ethernet                             |
| Device Dimension           | 9.6 x 17 x 10.5 inches                     |
| Device Weight              | 13 kg / 29 lbs.                            |
| Warranty                   | 1-year CARE Plus Plan included             |

# 3D PRINTING MADE EASY

# ONE

